

## **TOC Cleaning Validation Kit Simplifies Sampling for Cleaning Validation**

The TOC Cleaning Validation Kit from ITW Texwipe was designed to simplify sampling for cleaning validation. The TX3340 TOC Cleaning Validation kit contains the components necessary to sample 12 different areas in a production facility and provides for efficient transport of the samples from the production area to the laboratory. For facilities requiring a larger sampling kit, TX3342 TOC Cleaning Validation Bulk Kit is available with 72 vials and 144 swabs. The vials in the kit are compatible with the Anatel ANATOC™, Sievers Model 800, Tekmar Dohrman, OI Analytical Model 1010 and Shimadzu TOC-V Series Analyzers.

The kit offers several key advantages:

- TOC Cleaning Validation Kit provides a convenient method for transporting sampling vials and swabs to a sampling area
- Precleaned vials reduce TOC background and eliminate in-house processing costs
- Breakaway Large SnapSwab™ eliminates excessive handling and contamination

All kits offer the convenient packaging of the swabs and vials in a carrying box that enables efficient transport of the samples from production to laboratory with minimal chance of contamination. Moreover, the Large SnapSwab™ eliminates the need to cut the swab head." The internal head paddle provides firm support for sampling broad surfaces and flat areas. The breakaway handle allows the swab head to be placed in the vial with minimal handling.

The swabs and vials have been processed to reduce background TOC, with the Certificate of Analysis printed on the box. The double-bagged kit is also labeled with lot code and expiration date to simplify record keeping.